



intro to sample prep: SPE and SPME

Irina Galushko,
irina.galushko@merckgroup.com

Why Sample Prep?




May require a unique sample prep solution...



..but the same technology workflow for analysis




Time Spent on the Analytical Process



Sample prep

Sources of Chromatographic Errors



Source: LC/GC Magazin 2002



SPE


Solid Phase Extraction

MERCK


Real World & Real Samples

The Importance of Sample Preparation

Urine Sample without SPE



Urine Sample with SPE




Why is sample preparation required?

Collected Sample



GC, HPLC, or LC-MS/MS Analysis



Current Sample = Unsuitable for further analysis!!!... Why?

Too dirty- contains other sample matrix components that interfere with the analysis

Too dilute- analyte(s) not concentrated enough for quantitative detection

Present **sample matrix not compatible** with or harmful to the chromatographic column/system

SPE Formats


Sorbent particles held securely in place to withstand the force of the liquid flow.



96-well plates



Tubes



Disks



Loose/bulk sorbent (QuEChERS)



Online SPE



On the Inside: SPE Sorbents (Packing Materials)

The sorbent is the component of the tube responsible for the extraction. Most SPE sorbents are also used in HPLC applications, although with large particles in SPE. Some of the most common are:

Silica-based

- Reversed phase (C18, C8, cyano, phenyl)
- Normal phase (silica, diol, NH₂)
- Ion exchange (SAX, WCX, SCX)

Carbon-based

Polymer-based


- Various compositions
- Different functionalities

Others

- Florisil® (magnesium silicate)
- Alumina

Mixed-bed

- Combinations of nearly any of the above are possible in sequential layers



Supel™ Sphere dual-layer

SPE Strategies

There are 2 different elution strategies in SPE. Which one to choose depends on the goal of the extraction.


1. Bind-Elute Strategy

- Most common
- Bind: Analytes bind to tube, unwanted matrix components are washed off
- Elute: Eluant changed to remove analytes from tube
 - Different eluents can be used to fractionate the analytes
- Analytes are concentrated via evaporation prior to HPLC or GC analysis
- Sorbent types employing this: DSC-C18, Supel™-Select HLB, SupelMIP®, ENVI™-Carb Plus, PS/DVB, DSC-MCAX, ENVI™-Chrom P

2. Interference Removal Strategy

- Bind all unwanted matrix components and allow analytes to pass through during the sample loading stage
- Like chemical filtration
- Sorbent types employing this: HybridSPE®, QuEChERS, PSA, ENVI™-Carb, Dual Layer

Bind-Elute Strategy




Bind-Elute sorbent types:
DSC-C18, Supel™ -Select
HLB, SupelMIP®, ENVI™ -
Carb Plus, PS/DVB, DSC-
MCAX, ENVI™-Chrom P

Interference Removal Strategy

"Chemical Filtration"

Sample with Internal Standard in Matrix → Matrix adsorbed → Analytes & IS pass



Well-Established SPE Product Lines

Discovery®

Pharmaceutical focus

Tube and 96-well plates

Supelclean™ ENVI™

Environmental focus

ENVI™-Carb is a key product

Supel™-Select

Polymeric, “Universal SPE”

ENVI™-DSK™ disks

Porous glass fiber membranes embedded with sorbent particles



Sample Prep Star Products

Overcoated SPME

- Physically robust fiber for direct immersion that is less prone to chemical fouling.

Supelclean™ Ultra

- Dual layer cartridge for the cleanup of difficult matrices such as dry commodities (tea, spices, coffee, etc.)

Supel™MIP

- Molecularly imprinted polymers
- Highly selective for analytes in difficult matrices

Supel™ QuE

- QuEChERS tubes and supplies
- Pesticide Residue, PAH, PCB, PDBE analysis

Supel™ Tox

- Removes interferences associated with mycotoxin analysis

Supelclean™ EZ-POP NP

- Simple, effective extraction of lipophilic persistent organic pollutants (POPs) from oily samples

Supel™-Select Polymeric SPE

Supel-Select: What is Hydrophilic Polymer SPE?


Polymer chromatographic media designed for SPE

Comprises of a hydrophilic component and a hydrophobic component:

- Hydrophilic component examples:
 - N-vinyl pyrrolidone, methacrylate, hydroxyl, vinylamidizol
- Hydrophobic component examples:
 - Polystyrene, divinyl benzene


HLB Phase Chemistry:	Hydrophilic modified styrene polymer
SAX Phase Chemistry	Quaternary amine functionalized hydrophilic modified styrene polymer
SCX Phase Chemistry	Sulfonic acid functionalized hydrophilic modified styrene polymer
pH Compatibility:	0-14
Particle Size:	55-60 μm
MS Suitable:	Yes
Surface Area:	400-410 m^2/g
Pore Volume:	0.88 mL/g
Pore Size:	87 \AA

Unique Water-Wettable Oasis[®] HLB Copolymer




Supel-Select HLB: Statins from Rat Plasma

Supel-Select HLB SPE



Total Ion Chromatogram
(MRM, 4 pairs: 557.3/397.2) Rat
Plasma spiked with 5 ng/mL Statins



Absolute Recovery \pm RSD (n=3)

	5 ng/mL spike		100 ng/mL spike	
	Pravastatin	Atorvastatin	Pravastatin	Atorvastatin
Supel-Select HLB	84 \pm 8%	92 \pm 5%	103 \pm 4.2%	89 \pm 3.9%
Competitor W	83 \pm 17%	92 \pm 2%	104 \pm 2.2%	87 \pm 1.1%
Competitor P	77 \pm 5%	93 \pm 2%	102 \pm 3.0%	91 \pm 1.3%


Supel-Select HLB: Minimum Extractables

Assays today require greater sensitivity

SPE phase chemistry and hardware should impart minimum extractables

Each lot is tested for:

- Recovery
- LC-UV & LC-MS cleanliness
- Particle size
- Density
- Pore size
- Pore volume



Supel-Select Polymeric SPE

Why are they so popular in SPE?

- “Water – Wettable” – do not dry out => highly reproducible
- Amenable to generic methodology
- Often referred to as a universal SPE phase
- Can retain an extremely broad range of compounds (polar to non-polar; acidic – basic)
- Retained compounds easily eluted/desorbed with MeOH or similar solvent
- Reduces ion suppression in LC/MS
- Low UV and MS extractables
- 1000s of references using this technology

SUPELCO
Solutions within.™

Supel™-Select Polymeric SPE

Sample Prep Performance at the Perfect Price

Extract and Recover a Broad Range of Analytes from Aqueous Matrices

- Reduce Ion-Suppression
- Improve Reproducibility
- Enhance Sensitivity
- Save Money

SIGMA-ALDRICH

KZQ

www.sigmaaldrich.com/supel-select

MERCK

Current Supelclean™ Ultra 2400 Cartridge




- Cleanup **difficult matrices** prior to **pesticide residue analysis** by GC/MS/MS and LC/MS/MS
- **Dry commodities** (tea, spices, coffee, etc.), typically highly concentrated and with higher background than fresh samples
 - Pigments and oils
 - Not sufficiently cleaned by QuEChERS
- Dual layer SPE cartridge (1 mL and 3 mL) containing PSA, C18, specialty carbon, and Z-Sep
- **Specialized Carbon**
 - reduces pigmentation and allows for recovery of planar pesticides without toluene
- **Z-Sep sorbent**
 - Removes oils and some pigments, as was indicated in the cleanup of turmeric extracts for both GC and HPLC analysis

SupelMIP SPE – Molecularly Imprinted Polymer SPE

MIPs (molecularly imprinted polymers) are SPE products designed for the highly selective extraction of trace analytes from complex matrices

SupelMIP Phases and Methods Available for:


- PAHs in edible oil
- Non-steroidal anti-inflammatory drugs (NSAIDs) in wastewater and other sample matrices
- Nitroimidazoles in milk, eggs, and other food matrices
- Fluoroquinolones in bovine kidney, honey, and milk
- Chloramphenicol in milk, plasma, honey, urine and shrimp/prawns
- NNAL in urine
- TSNAs in urine and tobacco
- β -agonists in tissue, urine and wastewater
- Clenbuterol in urine
- Riboflavin in milk
- Patulin in fruit matrices
- Aminoglycosides in animal tissue, cell culture, and honey
- Bisphenol A from broth or milk-based matrices



- ❖ **Superior selectivity** => reduced ion-suppression => achieve lower detection limits
- ❖ **Robust & rapid methodology** => Save time, money, & headache
- ❖ **No method development req'd**

Supelclean™ EZ-POP NP:

For the extraction of non-polar compounds from edible oils



- Preferentially **retains fatty matrix interferences**
- **Non-polar compounds** are **washed through** and then analyzed
- **Compared to** other SPE methods:
 - Produces a **cleaner extract**
 - **Easier** and **more versatile methodology**
- The **final extracts** are in **acetonitrile**
 - Compatible with both **GC** and **HPLC**.



SPME

Solid Phase
Microextraction

Solid Phase Microextraction (SPME)

- Economical enrichment technique mainly for trace analysis
 - Semivolatiles & volatile (GC)
- Coated fused silica or metal fibers (adsorbent/particle & absorbent/film coatings)
- Initially for GC analysis, now extended to LC


Features:

- Very limited or no use of solvents
- All types of samples & matrixes
- Direct immersion or headspace
- Designs for manual, auto samplers and robots

Benefits:

- One-step extraction that is easy to automate
- Quantitative and reproducible Extractions
- Portable (field use) and reusable

[SigmaAldrich.com/spme](https://www.sigmaaldrich.com/spme)



Official Methods / Applications applying SPME Europe

ISO 27108 (former DIN 38407-34/Germany)

- Determination of **selected plant treatment agents and biocide products** - Method using solid-phase microextraction (SPME) followed by gas chromatography-mass spectrometry (GC-MS) – (Based on DIN 38407-F 34 from Germany)

ISO 17943 (former DIN 38407-41/Germany)

- German standard methods for the examination of water, waste water and sludge - Jointly determinable substances (group F) - Part 41: Determination of selected easily **volatile organic compounds in water** - Method using gas chromatography (GC-MS) after solid-phase micro extraction (SPME) (F 41)

OENORM A 1117, 2004-05-01

- Determination of **volatile compounds** in **cellulose-based materials** by Solid Phase Micro Extraction (SPME)

UNICIM 2237 / 2009 - Italian Method on SPME for air sampling

- Determinazione Delle Aldeidi Aerodisperse – Metodo per microestrazione su fase solida (SPME) ed analisi mediante gascromatografia accoppiata a spettrometria di massa (GC-MS)
- Determination of **air borne aldehydes by SPME/GC-MS** using **derivatisation on fiber** with PFBHA

Other official methods for SPME

US Methods

ASTM D 6438, 2005

- Standard Test Method for **Acetone, Methyl Acetate, and Parachlorobenzotrifluoride** Content of **Paints, and Coatings** by Solid Phase Microextraction-Gas Chromatography

ASTM D 6520, 2000

- Standard Practice for the Solid Phase Micro Extraction (SPME) of **Water** and its Headspace for the Analysis of **Volatile and Semi-Volatile Organic** Compounds

ASTM D 6889, 2003

- Standard Practice for Fast Screening for **Volatile Organic Compounds** in **Water** Using Solid Phase Microextraction (SPME)

ASTM E 2154, 2001

- Standard Practice for Separation and Concentration of **Ignitable Liquid Residues** from **Fire Debris** Samples by Passive Headspace Concentration with Solid Phase Microextraction (SPME)




EPA Method 8272 (Dec 2007)

- Parent and Alkyl **Polycyclic Aromatics** in **Sediment Pore Water** by SPME GC/MS

How SPME is used

1


Analyst has a **sample** they want to determine the content of



If sample prep is required


2

They **extract** the compounds from the sample onto the SPME fiber



3

They **desorb** the fiber into the GC instrument which tells them what was in the sample



GC chromatogram showing all the compounds that were extracted from the sample by the SPME fiber

The SPME Concept

SPME



Sample Adsorption

Please click on the numbered steps below for an animated sequence of the instruction.


- 1 Drill down septum piercing needle to avoid breakage
- 2 Insert needle into container
- 3 Adjust needle depth for aqueous sampling or headspace sampling
- 4 Extend plunger to expose fiber
- 5 Retract fiber before removing to avoid damaging the fiber.
- 6 Drill down septum piercing needle to avoid breakage.
- 7 Remove SPME Device

©2003, Sigma-Aldrich, All Rights Reserved

Sampling Technique: Headspace vs. Direct Immersion


Analytical considerations:

- Volatility of sample
- Extraction time concerns
- Sample matrix
- Selectivity of analytes



Adsorption Mechanism for SPME

SPME is quantitative!!!




SPME Fiber Coating: The Business End

An equilibrium is set up between analytes dissolved in the sample (solution or gas phase) and in the liquid coating on the fiber.

The fiber coating consists of:

- GC-type phases
- Particles



Types of SPME Fiber Coatings


Films – Absorption:

Coating	Type	Polarity
7 µm Polydimethylsiloxane (PDMS)	Absorbent	Nonpolar
30 µm PDMS	Absorbent	Nonpolar
100 µm PDMS	Absorbent	Nonpolar
85 µm Polyacrylate (PA)	Absorbent	Polar
60 µm PEG (Carbowax)	Absorbent	Polar

Particles – Adsorption:

Coating	Type	Polarity
85 µm Carboxen-PDMS	Adsorbent	Bipolar
65 µm PDMS-DVB	Adsorbent	Bipolar
55 µm/30 µm DVB/Carboxen-PDMS	Adsorbent	Bipolar

Molecular Weight Range for SPME Fibers



Effects of Salt and pH

- Salt usually increases analyte uptake
 - Use 25-30% NaCl to salt-out samples
 - Not necessary for large non-polar analytes, such as PAHs and large hydrocarbons, and may reduce recovery
- Lower pH to extract acidic compounds
- Raise pH to extract basic compounds
- Beware of stability of analytes at different pH levels

Thank You For Your Kind Attention!



Copyright © 2017 EMD Millipore Corporation. All Rights Reserved. MilliporeSigma, Millipore and the vibrant M are trademarks of Merck KGaA, Darmstadt, Germany. Supel, Supelclean, ENVI-Carb, and Stableflex are trademarks and Discovery, SLB, and SupelMIP are registered trademarks of Sigma-Aldrich Co. LLC. FocusLiner is a trademark of SGE Analytical Science Pty Ltd. Florisil is a registered trademark of U.S. Silica Company. Agilent is a registered trademark of Agilent Technologies, Inc.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.