



# intro to sample prep: SPE and SPME

Irina Galushko,  
[irina.galushko@merckgroup.com](mailto:irina.galushko@merckgroup.com)

# Why Sample Prep?




May require a unique sample prep solution...



..but the same technology workflow for analysis




# Time Spent on the Analytical Process



# Sample prep

## Sources of Chromatographic Errors



Source: LC/GC Magazin 2002



**SPE**


Solid Phase Extraction

**MERCK**


# Real World & Real Samples

## The Importance of Sample Preparation

Urine Sample without SPE



Urine Sample with SPE




# Why is sample preparation required?

Collected Sample



GC, HPLC, or LC-MS/MS Analysis



Current Sample = Unsuitable for further analysis!!!... Why?

**Too dirty**- contains other sample matrix components that interfere with the analysis

**Too dilute**- analyte(s) not concentrated enough for quantitative detection

Present **sample matrix not compatible** with or harmful to the chromatographic column/system

# SPE Formats


Sorbent particles held securely in place to withstand the force of the liquid flow.



96-well plates



Tubes



Disks



Loose/bulk sorbent (QuEChERS)



Online SPE





# On the Inside: SPE Sorbents (Packing Materials)

The sorbent is the component of the tube responsible for the extraction. Most SPE sorbents are also used in HPLC applications, although with large particles in SPE. Some of the most common are:

## Silica-based

- Reversed phase (C18, C8, cyano, phenyl)
- Normal phase (silica, diol, NH<sub>2</sub>)
- Ion exchange (SAX, WCX, SCX)

## Carbon-based

## Polymer-based


- Various compositions
- Different functionalities

## Others

- Florisil® (magnesium silicate)
- Alumina

## Mixed-bed

- Combinations of nearly any of the above are possible in sequential layers



Supel™ Sphere dual-layer

# SPE Strategies

**There are 2 different elution strategies in SPE. Which one to choose depends on the goal of the extraction.**


## 1. Bind-Elute Strategy

- Most common
- Bind: Analytes bind to tube, unwanted matrix components are washed off
- Elute: Eluant changed to remove analytes from tube
  - Different eluents can be used to fractionate the analytes
- Analytes are concentrated via evaporation prior to HPLC or GC analysis
- Sorbent types employing this: DSC-C18, Supel™-Select HLB, SupelMIP®, ENVI™-Carb Plus, PS/DVB, DSC-MCAX, ENVI™-Chrom P

## 2. Interference Removal Strategy

- Bind all unwanted matrix components and allow analytes to pass through during the sample loading stage
- Like chemical filtration
- Sorbent types employing this: HybridSPE®, QuEChERS, PSA, ENVI™-Carb, Dual Layer

# Bind-Elute Strategy




Bind-Elute sorbent types:  
DSC-C18, Supel™ -Select  
HLB, SupelMIP®, ENVI™ -  
Carb Plus, PS/DVB, DSC-  
MCAX, ENVI™-Chrom P

# Interference Removal Strategy

## "Chemical Filtration"

Sample with Internal Standard in Matrix → Matrix adsorbed → Analytes & IS pass



# Well-Established SPE Product Lines

## Discovery®

Pharmaceutical focus

Tube and 96-well plates

## Supelclean™ ENVI™

Environmental focus

ENVI™-Carb is a key product

## Supel™-Select

Polymeric, “Universal SPE”

## ENVI™-DSK™ disks

Porous glass fiber membranes embedded with sorbent particles





## Sample Prep Key Products for **Food Analysis**

### **Supelclean™ Ultra**

- Dual layer cartridge for the cleanup of difficult matrices such as dry commodities (tea, spices, coffee, etc.)

### **Supel™ MIP**

- Molecularly imprinted polymers
- Highly selective for analytes in difficult matrices

### **Supel™-Select Polymeric SPE**

- Yes, this can be used for food too!

### **Overcoated SPME**

- Physically robust fiber for direct immersion that is less prone to chemical fouling.

### **Supel™ QuE**

- QuEChERS tubes and supplies
- Pesticide Residue, PAH, PCB, PDBE analysis

### **Supel™ Tox**

- Removes interferences associated with mycotoxin analysis

### **Supelclean™ EZ-POP NP**

- Simple, effective extraction of lipophilic persistent organic pollutants (POPs) from oily samples

### **Supel™ Genie On-line SPE**

- Eliminates human error and reduces labor cost

# Supel-Select: What is Hydrophilic Polymer SPE?


## Polymer chromatographic media designed for SPE

### Comprises of a hydrophilic component and a hydrophobic component:

- Hydrophilic component examples:
  - N-vinyl pyrrolidone, methacrylate, hydroxyl, vinylamidizol
- Hydrophobic component examples:
  - Polystyrene, divinyl benzene


|                      |  |
|----------------------|--|
| HLB Phase Chemistry: | Hydrophilic modified styrene polymer                                 |
| SAX Phase Chemistry  | Quaternary amine functionalized hydrophilic modified styrene polymer |
| SCX Phase Chemistry  | Sulfonic acid functionalized hydrophilic modified styrene polymer    |
| pH Compatibility:    | 0-14   |
| Particle Size:       | 55-60 $\mu\text{m}$  |
| MS Suitable:         | Yes  |
| Surface Area:        | 400-410 $\text{m}^2/\text{g}$  |
| Pore Volume:         | 0.88 $\text{mL}/\text{g}$  |
| Pore Size:           | 87 $\text{\AA}$  |

Unique Water-Wettable Oasis<sup>®</sup> HLB Copolymer




# Supel-Select HLB: Statins from Rat Plasma

Supel-Select HLB SPE



Total Ion Chromatogram  
(MRM, 4 pairs: 557.3/397.2) Rat  
Plasma spiked with 5 ng/mL Statins



Absolute Recovery  $\pm$  RSD (n=3)

|                  | 5 ng/mL spike |              | 100 ng/mL spike |               |
|------------------|---------------|--------------|-----------------|---------------|
|                  | Pravastatin   | Atorvastatin | Pravastatin     | Atorvastatin  |
| Supel-Select HLB | 84 $\pm$ 8%   | 92 $\pm$ 5%  | 103 $\pm$ 4.2%  | 89 $\pm$ 3.9% |
| Competitor W     | 83 $\pm$ 17%  | 92 $\pm$ 2%  | 104 $\pm$ 2.2%  | 87 $\pm$ 1.1% |
| Competitor P     | 77 $\pm$ 5%   | 93 $\pm$ 2%  | 102 $\pm$ 3.0%  | 91 $\pm$ 1.3% |




# Supel-Select HLB: Minimum Extractables

Assays today require greater sensitivity

SPE phase chemistry and hardware should impart minimum extractables

Each lot is tested for:


- Recovery
- LC-UV & LC-MS cleanliness
- Particle size
- Density
- Pore size
- Pore volume



# Supel-Select Polymeric SPE

## Why are they so popular in SPE?

- “Water – Wettable” – do not dry out => highly reproducible
- Amenable to generic methodology
- Often referred to as a universal SPE phase
- Can retain an extremely broad range of compounds (polar to non-polar; acidic – basic)
- Retained compounds easily eluted/desorbed with MeOH or similar solvent
- Reduces ion suppression in LC/MS
- Low UV and MS extractables
- 1000s of references using this technology



The advertisement features a large, colorful, faceted sphere in the foreground, with a smaller, transparent, faceted sphere in the background. The background sphere is labeled 'SUPELCO' and 'SPE'. The text 'Supel™-Select Polymeric SPE' is prominently displayed in red, with the tagline 'Sample Prep Performance at the Perfect Price' below it. The Sigma-Aldrich logo is at the bottom right, and the text 'Extract and Recover a Broad Range of Analytes from Aqueous Matrices' is followed by a list of benefits: 'Reduce Ion-Suppression', 'Improve Reproducibility', 'Enhance Sensitivity', and 'Save Money'.

**SUPELCO**  
Solutions within.™

**Supel™-Select Polymeric SPE**  
Sample Prep Performance at the Perfect Price

Extract and Recover a Broad Range of Analytes from Aqueous Matrices  
Reduce Ion-Suppression  
Improve Reproducibility  
Enhance Sensitivity  
Save Money

**SIGMA-ALDRICH**

**KZQ**

[www.sigmaaldrich.com/supel-select](http://www.sigmaaldrich.com/supel-select)

**MERCK**

# Current Supelclean™ Ultra 2400 Cartridge



- Cleanup **difficult matrices** prior to **pesticide residue analysis** by GC/MS/MS and LC/MS/MS
- **Dry commodities** (tea, spices, coffee, etc.), typically highly concentrated and with higher background than fresh samples
  - Pigments and oils
  - Not sufficiently cleaned by QuEChERS
- Dual layer SPE cartridge (1 mL and 3 mL) containing:
  - **PSA** - removes acidic interferences
  - **C18** - retains some hydrophobic interferences
  - **Specialized Carbon** - reduces pigmentation and allows for recovery of planar pesticides without toluene
  - **Z-Sep sorbent** - removes oils and some pigments, as was indicated in the cleanup of turmeric extracts for both GC and HPLC analysis

[sigma-aldrich.com/supelcleanultra](http://sigma-aldrich.com/supelcleanultra)


MERCK

# SupelMIP SPE – Molecularly Imprinted Polymer SPE

MIPs (molecularly imprinted polymers) are SPE products designed for the highly selective extraction of trace analytes from complex matrices

## SupelMIP Phases and Methods Available for:


- PAHs in edible oil
- Non-steroidal anti-inflammatory drugs (NSAIDs) in wastewater and other sample matrices
- Nitroimidazoles in milk, eggs, and other food matrices
- Fluoroquinolones in bovine kidney, honey, and milk
- Chloramphenicol in milk, plasma, honey, urine and shrimp/prawns
- NNAL in urine
- TSNAs in urine and tobacco
- $\beta$ -agonists in tissue, urine and wastewater
- Clenbuterol in urine
- Riboflavin in milk
- Patulin in fruit matrices
- Aminoglycosides in animal tissue, cell culture, and honey
- Bisphenol A from broth or milk-based matrices



- ❖ Superior selectivity => reduced ion-suppression => achieve lower detection limits
- ❖ Robust & rapid methodology => Save time, money, & headache
- ❖ No method development req'd

# A New Approach: Analysis of Non-Polar POPs in Edible Oils

## Supelclean™ EZ-POP NP: Dual-layer SPE Cartridge Containing Florisil® and Z-Sep/C18 Mix



### Florisil® layer:

retains background constituents with polar functionality such as fatty acids




### Z-Sep/C18 layer:

retains fatty matrix through both Lewis acid/base and hydrophobic interactions

*\*Note: Yellow Color = Oil Matrix Removed From Sample*

- Easy sample preparation methodology using minimal volume of solvent.
- Final sample extracts compatible with GC or HPLC.

# Z-Sep/C18 And Proposed Lipid Retention Mechanism



# Supel™ Tox SPE for Mycotoxin Analysis

- Remove interferences associated with **mycotoxin analysis** in food and feed samples.
- **Quick, simple, and reproducible** sample cleanup solution
- **Compared to** the industry standard **Immunoaffinity Columns (IAC)**:
  - Decrease **sample prep time**
  - Increase **reproducibility**
  - More **user friendly**



Complemented by an extensive line of  
Certified Reference Materials (CRMs)

[sigma-aldrich.com/supeltox](http://sigma-aldrich.com/supeltox)  
[sigma-aldrich.com/mycotoxins](http://sigma-aldrich.com/mycotoxins)

- Supel Tox SPE Brochure (PFK)



**QUECHERS**

for Pesticide Residue  
Analysis

**MERCK**



# QuEChERS Methodology for Pesticide Residue Analysis....and beyond

- Stands for: **Q**uick, **E**asy, **C**heap **E**ffective, **R**ugged, **S**afe
- **Loose extraction salts** and **cleanup sorbents** in combination with **shaking** and **centrifugation** for purification
- Pesticide Residue, PAH, PCB, PDDBE analysis



## Official Methods:

- CEN Standard Method EN 15662
- AOAC Method 2007.01

|   |          |  |
|---|----------|--|
| 1 | SALTS    | Extraction: Phase separation and buffering<br>Cleanup: Water removal |
| 2 | SORBENTS | Cleanup: Remove unwanted interferences                               |
| 3 | TUBES    | Hold the ingredients and samples during extraction and clean up      |

[SigmaAldrich.com/quechers](http://SigmaAldrich.com/quechers)

# How Does the QuEChERS Method Work?

Weigh 10 g homogenized sample

Add 10 mL Acetonitrile (ACN)

Add Internal Standard

**Shake**

Add Salt Extraction Tube:  
AOAC 2007.01 = Acetate Tube  
EN15662:2008 = Citrate Tube

**Shake & centrifuge**

Transfer ACN layer to cleanup/sorbent tube

**Shake & centrifuge**

Adjust the pH if necessary

LC or GC Analysis



# QuEChERS: What's Inside – Cleanup tube

## SORBENTS

In AOAC and EN  
Official Methods

### PSA – “Primary Secondary Amine”

#### Supelclean™ PSA

- Removes: fatty acids, organic acids, polar pigments, sugars

### Carbon ENVI™-Carb

- Removes: chlorophyll, carotenoids

### C18 Discovery® DSC-18

- Removes: lipids, non-polar compounds

### Z-Sep Sorbents: Z-Sep, Z-Sep/C18, Z-Sep+

- Removes: lipids, pigments

Unique

Role of the sorbent is to trap interferences that co-extract with the analytes.

## QuEChERS Method: the choice of sorbent


| Interference | PSA      | C18      | C18/PSA  | ENVI-Carb | ENVI-Carb/ PSA | PSA/C18 /ENVI-Carb |
|--------------|----------|----------|----------|-----------|----------------|--------------------|
| Fats         |          | <b>X</b> | <b>X</b> |           |                | <b>X</b>           |
| Pigments     | <b>X</b> |          |          | <b>X</b>  | <b>X</b>       | <b>X</b>           |
| Sugars       | <b>X</b> |          | <b>X</b> |           | <b>X</b>       | <b>X</b>           |
| Acids        | <b>X</b> |          | <b>X</b> |           | <b>X</b>       | <b>X</b>           |

# Zirconia-based QuEChERS Sorbents Address Fatty Matrix Interferences





## Z-Sep Sorbents for QuEChERS

- Z-Sep/C18: For fatty and/or pigmented matrices containing < 15% fat
- Z-Sep+: For fatty and/or pigmented matrices containing > 15% fat
- Z-Sep: For the analysis of hydrophobic analytes in fatty matrices



# Proposed Retention Mechanism

## Hydrophobic Interactions



- Diminishes interferences from fatty matrices
- Removes various colors, specifically orange pigments
- Useful for analysis of pesticide residues
- Useful for analysis of PAH, PCB, PBDE and flame retardant analysis (specifically Z-Sep)

## Case Study:

### *Pesticides in Avocado by QuEChERS & GC/MS*

- Pesticide mix included hydrophobic compounds (e.g. organochlorines, hexachlorbenzene) and some other more polar classes **all GC/MS amenable**.



#### **Extraction and Cleanup Procedures**


1. Place 3 g of a homogenized avocado sample into a 50 mL centrifuge tube (Cat. No. 55248-U). Add spike solution if a spiked replicate.
2. Add 25 mL of acetonitrile and shake for one minute.
3. Add the contents of an Acetate Extraction Tube (Cat. No. 55234-U), and shake for one minute.
4. Centrifuge for five minutes.
5. Transfer 3 mL of the supernatant into the appropriate cleanup tube, Z-Sep+ (Cat. No. 55296-U) or PSA/C18 (Cat. No. 55229-U).
6. Shake for one minute, then centrifuge for three minutes.
7. Transfer 1 mL of the supernatant into an autosampler vial for GC/MS analysis.

## Case Study:

### *Pesticides in Avocado by QuEChERS & GC/MS*

Less color remained in  
the extract cleaned  
with  
Z-Sep+ vs. PSA/C18

#### GC/MS Background



#### Pigment Removal




GC-MS analysis of avocado  
extracts (scan mode) in the  
same y-scale



## Case Study:

### *Pesticides in Avocado by QuEChERS & GC/MS*

#### Results: Pesticide Recovery (Spiked samples each at 20 ng/g)



- Z-Sep+ showed better recovery overall.
- PSA/C18: matrix interference prevented analysis of cyfluthrin, cypermethrin and deltamethrin.
- Z-Sep+ showed better reproducibility than PSA/C18.



**SPME**

Solid Phase  
Microextraction

**MERCK**

# Solid Phase Microextraction (SPME)

- Economical enrichment technique mainly for trace analysis
  - Semivolatiles & volatile (GC)
- Coated fused silica or metal fibers (adsorbent/particle & absorbent/film coatings)
- Initially for GC analysis, now extended to LC


## Features:

- Very limited or no use of solvents
- All types of samples & matrixes
- Direct immersion or headspace
- Designs for manual, auto samplers and robots

## Benefits:

- One-step extraction that is easy to automate
- Quantitative and reproducible Extractions
- Portable (field use) and reusable


[SigmaAldrich.com/spme](https://www.sigmaaldrich.com/spme)



# How SPME is used

1


Analyst has a **sample** they want to determine the content of



If sample prep is required


2

They **extract** the compounds from the sample onto the SPME fiber



3

They **desorb** the fiber into the GC instrument which tells them what was in the sample



GC chromatogram showing all the compounds that were extracted from the sample by the SPME fiber

# The SPME Concept

**SPME**



## Sample Adsorption

Please click on the numbered steps below for an animated sequence of the instruction.


- 1 Drill down septum piercing needle to avoid breakage
- 2 Insert needle into container
- 3 Adjust needle depth for aqueous sampling or headspace sampling
- 4 Extend plunger to expose fiber
- 5 Retract fiber before removing to avoid damaging the fiber.
- 6 Drill down septum piercing needle to avoid breakage.
- 7 Remove SPME Device

©2003, Sigma-Aldrich, All Rights Reserved

# Sampling Technique: Headspace vs. Direct Immersion


Analytical considerations:

- Volatility of sample
- Sufficient vapor pressure
- Extraction time concerns
- Sample matrix
- Selectivity of analytes



# Adsorption Mechanism for SPME

**SPME is quantitative!!!**




# SPME Fiber Coating: The Business End

An equilibrium is set up between analytes dissolved in the sample (solution or gas phase) and in the liquid coating on the fiber.

The fiber coating consists of:

- GC-type phases
- Particles





# Types of SPME Fiber Coatings


## Films – Absorption:

| Coating                          | Type      | Polarity |
|----------------------------------|-----------|----------|
| 7 µm Polydimethylsiloxane (PDMS) | Absorbent | Nonpolar |
| 30 µm PDMS                       | Absorbent | Nonpolar |
| 100 µm PDMS                      | Absorbent | Nonpolar |
| 85 µm Polyacrylate (PA)          | Absorbent | Polar    |
| 60 µm PEG (Carbowax)             | Absorbent | Polar    |

## Particles – Adsorption:

| Coating                       | Type      | Polarity |
|-------------------------------|-----------|----------|
| 85 µm Carboxen-PDMS           | Adsorbent | Bipolar  |
| 65 µm PDMS-DVB                | Adsorbent | Bipolar  |
| 55 µm/30 µm DVB/Carboxen-PDMS | Adsorbent | Bipolar  |

# Molecular Weight Range for SPME Fibers



# Effects of Salt and pH

- Salt usually increases analyte uptake
  - Use 25-30% NaCl to salt-out samples
  - Not necessary for large non-polar analytes, such as PAHs and large hydrocarbons, and may reduce recovery
- Lower pH to extract acidic compounds
- Raise pH to extract basic compounds
- Beware of stability of analytes at different pH levels

# Thank You For Your Kind Attention!



Copyright © 2017 EMD Millipore Corporation. All Rights Reserved. MilliporeSigma, Millipore and the vibrant M are trademarks of Merck KGaA, Darmstadt, Germany. Supel, Supelclean, ENVI-Carb, and Stableflex are trademarks and Discovery, SLB, and SupelMIP are registered trademarks of Sigma-Aldrich Co. LLC. FocusLiner is a trademark of SGE Analytical Science Pty Ltd. Florisil is a registered trademark of U.S. Silica Company. Agilent is a registered trademark of Agilent Technologies, Inc.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.