

---

# Supel™ QuE

Multi-Residue Pesticide Analysis  
in Food and Agricultural Products  
using the QuEChERS Method.

---



# SUPEL™ QuE

## Multi-residue Pesticide Analysis in Food and Agricultural Products

In “QuEChERS” methodology, the use of loose extraction salts and cleanup sorbents in combination with shaking and centrifugation results in a **Quick, Easy, Cheap, Effective, Rugged and Safe** sample cleanup technique. The “QuEChERS” method has emerged as a sample prep technique popular in the area of multi-residue pesticide analysis in food and agricultural products, and is formalized in method EN15662:2008 and AOAC 2007.01.<sup>1-2</sup> Recently, QuEChERS has been expanded into areas including PAH, PCB, PBDE and flame retardant analysis.<sup>3</sup>

The Supel QuE line of centrifuge tubes contains predetermined amounts of salts and SPE sorbents to support the most common method configurations used today for QuEChERS.

### Features and Impact

QuEChERS significantly improves laboratory efficiency and throughput. This procedure requires only a few quantities of solvent and is capable of generating recoveries of 70-120% with RSDs <5% for a wide range of compounds.



## A PROCESS IN TWO STEPS

### 1-Extraction

- Solvent extraction techniques are designed to yield as much analytes from the base matrix as possible
- Solvent selection is important to minimize co-extracting compounds
- Analytes are extracted from the matrix with acetonitrile and salts/buffers

### 2-Cleanup

- Sample cleanup is necessary to reduce interferences
- Interferences can damage analytical instrumentation and complicate analyte identification and quantification
- The use of PSA, C18, GCB (Graphitized Carbon Black) or Z-Sep allow removal of sugars, lipids, sterols, organic acids, proteins, carotenoids, chlorophyll and other pigments prior to a GCMS-MS or and LCMS-MS analysis

## OFFICIAL QUECHERS METHODS

### EN 15662:2008

- The European method includes sodium chloride to limit polar interferences and several buffering citrate based reagents to preserve base sensitive analytes
- The use of sodium hydroxide in the citrate step should be avoided as it damages the sorbent used in the cleanup step

### AOAC 2007.01

- Employs 1% acetic acid in acetonitrile and sodium acetate buffer to protect base sensitive analytes from degradation
- A USDA study has demonstrated that this method provides superior recovery for pH sensitive compounds when compared to the other QuEChERS methods. Since the approach uses acetic acid in the extraction step, the PSA sorbent used in the cleanup step might get overloaded making it less effective in cleanup and possibly causing GC resolution/background issues.

#### EN 15662:2008

#### 2007.01 AOAC

#### Step 1 - Sample Preparation

Freeze samples to -20 °C  
Homogenize with dry ice until a free flowing powder is formed

Freeze samples to -20 °C  
Homogenize with dry ice until a free flowing powder is formed

#### Step 2 - Sample Extraction

Weight 10 g sample into a 50 mL centrifuge tube

Transfer 10-15 g homogenized food sample to 50 mL PTFE centrifuge tube

Add 10 mL ACN + 100 µL I.S. solution. Shake vigorously for 1 minute

For each 15 g sample, add 15 mL 1% acetic acid in ACN + contents of Supel QuE Acetate (AC) Tube (55234-U) + 75µL I.S. solution.

Add contents of Supel™ QuE Citrate Extraction Tube (55227-U).  
Shake for 1 minute and centrifuge for 5 minute at 3000 U/min

Shake for 1 minute and centrifuge for 1 minute at 1500 U/min

#### Step 3 - Cleanup - QuEChERS Technique

Transfer an aliquot of the ACN layer (supernatant) to a dispersive cleanup tube. Shake for 30 seconds (2 minute if ENVI-Carb™)  
Centrifuge 5 minute at 3000 U/min

Transfer an aliquot of the ACN layer (supernatant) to a dispersive cleanup tube. Shake for 30 seconds (2 minute if ENVI-Carb)  
Centrifuge 1 minute at 1500 U/min

Acidify an aliquot of supernatant with formic acid 5% in ACN  
(10 µL to every mL of supernatant isolated)

Further processing may be necessary prior to chromatographic analysis  
(addition of formic acid, evaporation and reconstitution with toluene)

## QUECHERS: BEYOND THE BASICS

QuEChERS was originally designed for fruits and vegetables and proven to work well especially for high water content fruits and vegetables (>80%) under slightly acidic conditions (pH 5-6).<sup>1-4</sup> However, there are commodities that do not fit into this food type. Adjustments to the typical QuEChERS procedures extend its usage outside of these typical or easy commodities.

To improve the extraction efficiency of low moisture containing commodities, addition of water to the samples prior to extraction is recommended (see table to the right).

To reduce co-extracted materials from acidic foods,<sup>5</sup> buffering is needed. The buffering of the EN method is intended to produce this pH 5-6 condition during extraction.<sup>3</sup> The buffering capacity of existing methods cannot adequately correct for the low pH of citrus fruits (with pH<3). In this case the EN method calls for addition of 5 N (or M, mol/L) sodium hydroxide solution to the salt mixes. In case of lemons and lime 600 µL are recommended.<sup>3</sup> For commodities like raspberry, addition of 200 µL to the extraction tube is recommended. Acidic foods outside the citrus family can be pH adjusted and generally produce comparable high-quality data as expected with non acidic foods.<sup>1-4</sup>

Sample Types	Sample Weight	Water Added
Fruits and vegetables > 80% water content	10 g	—
Fruits and vegetables 25-80% water content	10 g	x g
Cereals	5 g	10 g
Dried fruits	5 g	7,5 g
Honey	5 g	10 g
Spices	2 g	10 g

x = 10g – water amount in 10 g sample.


### References

1. Fast and Easy Multiresidue Method Employing Acetonitrile Extraction/Partitioning and "Dispersive Solid-Phase Extraction" for the Determination of Pesticide Residues in Produce. (M. Anastassiades, S.J. Lehotay, D. Stajnbauer, F.J. Schenck, J. AOAC International 86 (2003) 412).
2. Pesticide Residues in Foods by Acetonitrile Extraction and Partitioning with Magnesium Sulfate. (AOAC Official Method 2007.01).
3. Foods of Plant Origin—Determination of Pesticide Residues Using GC-MS and/or LC-MS/MS Following Acetonitrile Extraction/Partitioning and Clean-up by Dispersive SPE (QuEChERS-method). (EN 15662 Version 2008).
4. Comparison of QuEChERS sample preparation methods for the analysis of pesticide residues in fruits and vegetables. (S. Lehotay, K. Son, H. Kwon, U. Koesukwiwat, W. Fu, K. Mastovska, E. Hoh, N. Leepipatpiboon, J. Chromatogr. A, 1217 (2010) 2548).
5. The QuEChERS Method –Background Information and Recent Developments, Community Reference Laboratory Pesticide Residues using Single Residue Methods. (M. Anastassiades, 1st Joint CRL-Workshop –Stuttgart, presentation (2006)).

### Supel™ QuE Verde for Challenging Compounds in Green Matrices – NEW!

Supel QuE Verde for QuEChERS combines a novel carbon with zirconia coated silica (Z-Sep+) to provide an optimum balance between analyte recovery and color removal. This sorbent combination has been shown to provide recoveries in the range of 70% to 120% of even the most challenging planar pesticides while maintaining >95% pigment removal in high chlorophyll matrices.


Supel QuE Verde is a mixture of an improved graphitized carbon black (GCB), Z-Sep+, and primary-secondary amine (PSA). The improved GCB has been optimized to balance chlorophyll removal and improve recoveries of planar pesticides. Z-Sep+ is a silica that is functionalized with both zirconia and C18. Zirconia will retain some fats and carotenoids, while C18 retains hydrophobic interferences. The PSA in the mix functions to remove acidic interferences. When used to clean samples containing chlorophyll, this sorbent blend will provide better recovery of planar pesticides than sorbents containing traditional GCB.



Average Percent Recoveries From Spinach and Oregano Extracts Spiked at 50 ng/mL After Cleanup With the 2 mL Supel QuE Verde Tube (n=3).

Discover Supel QuE Verde on [sigma-aldrich.com/verde](http://sigma-aldrich.com/verde)

## Classification of Fruits and Vegetables (as defined in the CEN Guidelines)

Commodity Groups	Typical Commodity Categories	Typical Representative Commodities
<b>1. High water content</b> 	Pome fruit Stone fruit Other fruit Alliums Fruiting vegetables/cucurbits Brassica vegetables Leafy vegetables and fresh herbs Stem and stalk vegetables Forage/fodder crops Fresh legume vegetables Leaves of root and tuber vegetables Fresh fungi Root and tuber vegetables or feed	Apples, pears Apricots, cherries, peaches, Bananas Onions, leeks Tomatoes, peppers, cucumber, melon Cauliflower, Brussels-sprouts, cabbage, broccoli Lettuce, spinach, basil Celery, asparagus Fresh alfalfa, fodder vetch, fresh sugar beets Fresh peas with pods, peas, mange tout, broad beans, runner beans, French beans Sugar beet and fodder beet tops Champignons, canterelles Sugar beet and fodder beet roots, carrots, potatoes, sweet potatoes
<b>2. High acid content and high water content<sup>1</sup></b> 	Citrus fruit Small fruit and berries Other	Lemons, mandarins, tangerines, oranges Strawberry, blueberry, raspberry, black currant, red currant, white currant, grapes Kiwifruit, pineapple, rhubarb
<b>3. High sugar and low water content<sup>2</sup></b> 	Honey, dried fruit	Honey, raisins, dried apricots, dried plums, fruit jams
<b>4a. High oil content and very low water content</b> 	Tree nuts Oil seeds Pastes of tree nuts and oil seeds Oils from tree nuts, oil seeds and oily fruits	Walnuts, hazelnuts Oilseed rape, sunflower, cotton-seed, soybeans, peanuts, sesame, etc. Peanut butter, tahina, hazelnut paste Olive oil, rapeseed oil, sunflower oil, pumpkin seed oil
<b>4b. High oil content and intermediate water content</b> 	Oily fruits and products	Olives, avocados and pastes thereof
<b>5. High starch and/or protein content and low water and fat content</b> 	Dry legume vegetables/pulses Cereal grain and products thereof	Field bean, dried broad bean, dried haricot bean (yellow, white/navy, brown, speckled), lentils Wheat, rye, barley and oat grain; maize, rice, wholemeal bread, white bread, crackers, breakfast cereals, pasta
<b>6. Difficult or unique commodities</b> 	—	Hops, cocoa beans and products thereof, coffee, tea spices

<sup>1</sup> If a buffer is used to stabilize the pH changes in the extraction step, then commodity Group 2 can be merged with commodity Group 1.

<sup>2</sup> Where commodities of Group 3 are mixed with water prior to extraction to achieve a water content of >70%, this commodity group may be merged with Group 1. The RLs should be adjusted to account for smaller sample portions (e.g. if 10 g portions are used for commodities of Group 1 and 5 g for Group 3, the RL of Group 3 should be twice the RL of Group 1 unless a commodity belonging to Group 3 is successfully validated at a lower level).

# HOW TO CHOOSE A SPE SORBENT?

Supelco offers centrifuge tubes in bulk quantities, each packed with predetermined amounts of salts and SPE sorbents to support the most common method configurations used today. This saves the analyst the time taken to weigh out the individual materials multiple times.

## Cleanup Products and Their Characteristics

### Zirconia on Silica (Z, Zep, Z-Sep/C18, Z-Sep+)

- Lewis Acid-Lewis base interactions between hydroxyl groups and Zirconium
- Hydrophobic interactions between triglycerides and C18 (Z-Sep/C18, Z-Sep+)
- Significantly diminishes fatty matrix interferences and various pigments
- Provides more robust LC/MS and GC/MS methods by eliminating problematic matrix interferences
- Can replace C18 and PSA phases in current methods without additional method development

### Primary Secondary Amine (PSA)

- Polymerically bonded, ethylenediamine-N-propyl phase that contains both primary and secondary amines
- A weak anion exchanger with a pKa of 10.1 and 10.9
- Similar to aminopropyl SPE phases (NH<sub>2</sub>) in terms of selectivity, but has a much higher capacity due to presence of secondary amine (0.98-1.05 meq/g)
- Has been shown to significantly reduce matrix-enhancement effects encountered during the GC analysis of food products
- Bidentate nature of ligands allow for chelation

### ENVI-Carb

- Extreme affinity for organic polar and nonpolar compounds with ringstructures/conjugated double bonds from both nonpolar and polar matrices when used under reversed-phase conditions
- Carbon surface comprised of hexagonal ring structures, interconnected and layered into graphitic sheets
- Non-porous; nature of the carbon phase allows for rapid processing; adsorption does not require analyte diffusion into pores


Zirconia
<ul style="list-style-type: none"><li>• Removes pigments and lipids or fats (content below 15% fats use Z-Sep/C18, for above 15% use Z-Sep+, for hydrophobic analytes use Z-Sep)</li></ul>

PSA
<ul style="list-style-type: none"><li>• Used in the removal of sugars and fatty acids, organic acids, lipids and polar pigments</li><li>• When used in combination with C18, additional lipids and sterols can be removed</li><li>• When used in combination with GCBs, higher amounts of pigments can be removed</li></ul>

ENVI-Carb™
<ul style="list-style-type: none"><li>• Strong sorbent for removing pigments, polyphenols, and other polar compounds</li></ul>

C18
<ul style="list-style-type: none"><li>• Removes long carbon chain/fatty compounds, sterols and other nonpolar interferences</li></ul>

### Retention Mechanism for Fats on Supel QuE Z-Sep Sorbents



## WHY CHOOSE SUPEL QUE?

### We were among the first to introduce QuEChERS to the market

- Excellent knowledge of the technique and manufacturing
- Extended choice of adsorbents and tubes sizes (incl. 15 mL for automated shaker)

### Innovative new adsorbents with unique properties

- Z-Sep sorbents for fatty and difficult food
- Supel QuE Verde for planar pesticides in green matrices

### MgSO<sub>4</sub>

- Supelco heats the MgSO<sub>4</sub> to 550 °C prior to packing to ensure dryness

### Customization Services

- Adsorbent mixes and configuration for special requirements
- Hardware options (PP containers, glass vials)
- Components to be packed (salts, adsorbents)



### Did you know . . .

The Sigma-Aldrich® pesticide standards product line is the most comprehensive portfolio available in the market. Under the well-known brands PESTANAL® and TraceCERT®, we proudly offer more than 1300 high purity pesticide and pesticide metabolite standards and certified reference materials for food and environmental analysis including:

- Certified Reference Materials (CRMs) TraceCERT
- CRM Matrix Standards
- Isotope Labeled Standards
- Pesticide Metabolite Standards

We add new products to the portfolio regularly, in order to keep up-to-date with new developments in pesticide analysis.

For more information, visit  
[sigma-aldrich.com/pesticides](http://sigma-aldrich.com/pesticides)



*Coming soon: Pesticide residue analysis iPad App, featuring sampling and analysis methods. Watch the pesticide website for the announcement.*

# SUPEL QUE PRODUCTS FOR QUECHERS AND RELATED PRODUCTS

## Products for the Extraction Step

### Supel QuE Tubes for the Extraction Step

Description	Qty	Cat. No.
<b>Method EN15662 :2008 (12 mL centrifuge tubes)</b>		
Supel QuE Citrate (EN) Tube – 4 g MgSO <sub>4</sub> , 1 g NaCl, 0.5 g NaCitrate dibasic sesquihydrate, 1 g NaCitrate tribasic dihydrate	50	55227-U
Supel QuE Citrate/Sodium Bicarbonate (EN) Tube – 4 g MgSO <sub>4</sub> , 5 g NaBicarbonate, 1 g NaCl, 0.5 g NaCitrate dibasic	50	55237-U
<b>Method AOAC 2007.01 (12 mL centrifuge tubes)</b>		
Supel QuE Acetate (AC) Tube – 6 g MgSO <sub>4</sub> , 1.5 g NaAcetate	50	55234-U
<b>Non buffered extraction tubes (12 mL centrifuge tubes)</b>		
Supel QuE Non-Buffered Tube 1 – 4 g MgSO <sub>4</sub> , 1 g NaCl	50	55294-U
Supel QuE Non-Buffered Tube 2 – 6 g MgSO <sub>4</sub> , 1.5 g NaCl	50	55295-U

### QuEChERS Shakers and Accessories

Description	Qty	Cat. No.
<b>Benchmark Benchmixer™ XL Laboratory Shakers</b>		
QuEChERS Shaker and Rack Starter Kit, USA compatible plug, AC input 115 V	—	55278-U
QuEChERS Shaker and Rack Starter Kit, EU compatible, schuko plug, AC input 230 V	—	55438-U
Multi-tube Vortexer, USA compatible plug, AC input 115 V	1 ea.	Z765503
Multi-tube Vortexer, EU compatible Schuko plug, AC input 230 V	—	Z765511
<b>Benchmark Benchmixer XL Laboratory Shaker Racks</b>		
50 mL QuEChERS Extraction Tube Shaker Rack	pk 1	55279-U
15 mL QuEChERS Cleanup Tube Shaker Rack	pk 1	Z765589
2 mL QuEChERS Cleanup Tube Shaker Rack	pk 1	Z765554

### Related Products

Description	Qty	Cat. No.
<b>Ascentis® Express HPLC Columns (2.7 µm particles)</b>		
RP-Amide, 10 cm × 2.1 mm I.D.	1	53913-U
<b>Capillary GC Columns</b>		
SLB-5ms, 30 m × 0.25 mm I.D. × 0.25 µm	1	28471-U

You want to know more about our columns or you are interested in other sampling techniques? Request our technical brochures

- GC column selection guide (KCX)
- SPE brochure (FEB)
- Ascentis Express HPLC columns (JHD)





## Products for Cleanup (QuEChERS) Step

### Method – EN15662:2008

Matrices*	Product Description	Cat. No. 2 mL Tubes	Cat. No. 15 mL Tubes	
Group 1, 2 and 3	General fruits and vegetables <i>lightly pigmented</i>	150 mg Supelclean™ PSA, 900 mg MgSO <sub>4</sub>	—	55437-U
		25 mg Supelclean PSA, 150 mg MgSO <sub>4</sub>	55172-U	—
	Pigmented fruits and Vegetables <i>Moderate levels of chlorophyll and carotinoides</i>	150 mg Supelclean PSA, 15 mg Supelclean ENVI-Carb, 900 mg MgSO <sub>4</sub>	—	55446-U
		25 mg Supelclean PSA, 2,5 mg Supelclean ENVI-Carb, 150 mg MgSO <sub>4</sub>	55147-U	—
	Pigmented fruits and vegetables <i>High levels of chlorophyll and carotinoides</i>	150 mg Supelclean PSA, 45 mg Supelclean ENVI-Carb, 900 mg MgSO <sub>4</sub>	—	55464-U
		25 mg Supelclean PSA, 7,5 mg Supelclean ENVI-Carb, 150 mg MgSO <sub>4</sub>	55176-U	—
Group 4, 5 and 6	Pigmented fruits and vegetables with waxes/lipids High Lipid Content	150 mg Supelclean PSA, 150 mg Discovery® DSC-18, 900 mg MgSO <sub>4</sub>	—	55439-U
		25 mg Supelclean PSA, 25 mg Discovery DSC-18, 150 mg MgSO <sub>4</sub>	55173-U	—

### Method – AOAC 2007.01

Matrices*	Product Description	Cat. No. 2 mL Tubes	Cat. No. 15 mL Tubes	
Group 1, 2 and 3	General fruits and vegetables <i>lightly pigmented</i>	400 mg Supelclean PSA, 1200 mg MgSO <sub>4</sub>	—	55466-U
		50 mg Supelclean PSA, 150 mg MgSO <sub>4</sub>	55287-U	—
	Pigmented fruits and vegetables <i>Moderate levels of chlorophyll and carotinoides</i>	50 mg Supelclean PSA, 150 mg MgSO <sub>4</sub> , 50 mg ENVI-Carb	on request	—
		50 mg Supelclean PSA, 150 mg MgSO <sub>4</sub> , 50 mg Discovery DSC-18, 50 mg ENVI-Carb	55289-U	—
	Pigmented fruits and Vegetables <i>High levels of chlorophyll and carotinoides</i>	400 mg Supelclean PSA, 1200 mg MgSO <sub>4</sub> , 400 mg Discovery DSC-18, 400 mg ENVI-Carb	—	55474-U
		50 mg Supelclean PSA, 150 mg MgSO <sub>4</sub> , 50 mg Discovery DSC-18	55288-U	—
Group 4, 5 and 6	Pigmented fruits and vegetables with waxes/lipids High Lipid Content	400 mg Supelclean PSA, 1200 mg MgSO <sub>4</sub> , 400 mg Discovery DSC-18	—	55470-U

\* Please contact your local Technical Service Team

### Specialty Products for Method EN15662 or AOAC 2007.01 – Alternatives

Matrices*	Colonne1	Product description	Cat.No. 2 mL tubes	Cat.No. 15 mL tubes
Group 4, 5 and 6	Fatty or pigmented matrix <i>With less than 15% fat</i>	20 mg Z-Sep, 50 mg Discovery DSC-C18 (alternative)	55284-U	—
		120 mg Z-Sep, 300 mg C18 (alternative)	—	55506-U
	<i>Hydrophobic analytes in fatty matrices</i>	75 mg Z-Sep (alternative)	55411-U	—
		50 mg Z-Sep, 150 mg MgSO <sub>4</sub> (alternative)	55417-U	—
		500 mg Z-Sep, (alternative)	—	55491-U
		300 mg Z-Sep, 900 mg MgSO <sub>4</sub> (alternative)	—	55503-U
	Fatty matrices <i>With greater than 15% fat</i>	75 mg Z-Sep+ (alternative)	55408-U	—
		50 mg Z-Sep+, 150 mg MgSO <sub>4</sub> (alternative)	55414-U	—
		500 mg Z-Sep+, (alternative)	—	55486-U
		300 mg Z-Sep+, 900 mg MgSO <sub>4</sub> (alternative)	—	55511-U
<b>NEW</b>	Improved recovery of <i>planar pesticides</i> in green matrices	Supel QuE Verde, 150 mg MgSO <sub>4</sub> , 10 mg ENVI-Carb Y, 50 mg PSA, 60 mg Z-Sep+	55447-U	—
		Supel QuE Verde, 1200 mg MgSO <sub>4</sub> , 80 mg ENVI-Carb Y, 400 mg PSA, 480 mg Z-Sep+	—	55442-U

\*See page 5 for group definitions.



## APPLICATIONS

### Intralab Validation of the EN 15662 Method for the Determination of 200 Pesticide Residues Using a Fused-Core® Ascentis® Express RP-Amide HPLC Column by LC-MS/MS and Cleanup by Dispersive SPE (QuEChERS)

Enio Belotti, Luca Meni, Marco Ruggeri, Water and Life Entratico (BG) Italy

## Experimental

The European guideline SANCO/3131/2007 document (Method Validation and Quality Control Procedures for Pesticide Residues Analysis in Foods and Feeds) was followed for the following different representative fruit and plant origin matrices:

- Pear as sugar matrix
- Kiwi as acid matrix
- Lettuce as chlorophyll matrix
- Maize flour as cereal

The QuEChERS Multi-residue method described in EN 15662 was used to prepare the final samples to be injected into the LC/MS/MS system.

## Chromatography

The HPLC method developed here utilizes an Ascentis Express RP-Amide (Supelco®, Bellefonte) as an alternative to the C18 column traditionally used in the EN 15662 method.

Combining an embedded polar group (EPG) stationary phase with the Fused-Core particles, Ascentis Express RP-Amide provides a host of useful benefits to the HPLC chromatographer.

The benefits come from both the phase technology and the particle technology, and can be summarized as:

### Fused-Core Benefits

- Twice the efficiency of traditional 3 µm HPLC columns
- Half the backpressure of sub-2 micron columns
- Capable of UHPLC performance on traditional HPLC systems

### RP-Amide Benefits

- Alternative reversed-phase selectivity to C18
- Improved peak shape for bases
- 100% aqueous compatible reversed-phase column

### Alternative Selectivity


Ascentis Express RP-Amide provides increased selectivity for polar compounds, especially those that can act as a hydrogen-bond donor. The selectivity differences between the RP-Amide and the C18 can be a useful tool in method development. In many cases, when peaks co-elute on a C18 phase, the RP-Amide can be substituted to achieve separation without a change in mobile phase.

For more details, download the Reporter Europe Vol. 39 /2010 at [sigma-aldrich.com/thereporter](http://sigma-aldrich.com/thereporter)

### Operating conditions:

column: Ascentis Express RP Amide,  
10 cm x 2.1 mm ID, 2.7 µm  
particle size  
HPLC: Shimadzu  
Prominence UFLC XR  
MS/MS detector: Applied Biosystems  
API 3200  
mobile phase A: Ammonium formate  
solution in water  
(5 mmol/l, 0.1% formic acid)  
mobile phase B: Ammonium formate  
solution in methanol  
(5 mmol/l, 0.1% formic acid)  
column temp: 40 °C  
injection volume: 2 µL  
elutin: gradient

Time	Mobile Phase A %	Mobile Phase B %
0	95	5
0.5	90	10
12	5	95
15	85	95





# Determination of 113 pesticide residues in high oil vegetal commodities

Lukasz Rajski, Ana Lozano, Ana Uclés, Carmen Ferrer, Amadeo R. Fernández-Alba - *J. Chromatogr. A* 1304 (2013) 109– 120

Prof. Fernandez-Alba\* and his team have evaluated several extraction methods in term of recoveries and extraction precision for 113 pesticides in avocado including QuEChERS with various dSPE cleanups, miniLuke, and ethyl acetate. This work is described in a paper published in *Journal of Chromatography A*.

Supel QuE Z-Sep and Z-Sep+, two sorbents containing  $ZrO_2$ , were used to improve fat removal from the extracts. The findings were that the QuEChERS protocol with Z-Sep cleanup showed the highest number of pesticides with recoveries in the 70–120% range, along with the lowest amount of co-extracted matrix compounds for avocado and almond matrices. As part of method validation, recoveries at two levels (10 and 50  $\mu\text{g}/\text{kg}$ ), limit of quantitation, linearity, matrix effects, as well as the inter- and intra-day precision were studied.

In addition to providing the best recoveries, Z-Sep cleaned extracts also exhibited the lowest average % RSD values.


The authors also showed that apart from the application of a suitable sorbent, the method for sample reconstitution after evaporation (if applied) is very important. Reconstitution in a mixture of 1:9 (v/v) acetonitrile: water was inefficient. In avocado, in order to achieve good recoveries for a large number of pesticides, it was necessary to increase the amount of acetonitrile to 30%. Almond samples, on the other hand, required 40% acetonitrile.

Even with these percentages of acetonitrile, sediment was observed in the extracts cleaned with PSA/C18, while the Z-Sep cleaned extracts were clear. Almonds were the more difficult matrix to analyze. Out of 113 pesticides analyzed, at a spiking level of 50  $\mu\text{g}/\text{kg}$ , 94 had recoveries ranging from 70 to 120%. At a spiking level of 10  $\mu\text{g}/\text{kg}$ , 92 pesticides had recoveries in this range. From avocado, 107 of the 113 pesticides had recoveries of 70–120% at both 10 and 50  $\mu\text{g}/\text{kg}$  spiking levels.


\* Pesticide Residue Research Group, European Union Reference Laboratory (EURL), Department of Hydrogeology and Analytical Chemistry, University of Almería, Ctra. Sacramento S/N, La Cañada de San Urbano, 04120, Almería, Spain

To find the abstract & link to this article, see references at [sigma-aldrich.com/zsep](http://sigma-aldrich.com/zsep)

(a) Percentage of total number of evaluated pesticides with recoveries in the range of 70 to 120%, in avocado. Red bars: evaporation of 0.50 mL, reconstitution in 0.05 mL of ACN + 0.45 mL of H<sub>2</sub>O; blue bars: evaporation of 0.50 mL, reconstitution in 0.30 mL of ACN + 0.70 mL of H<sub>2</sub>O.



(b) Percentage of total number of evaluated pesticides with recoveries in the range of 70 to 120% in almonds using Z-Sep sorbent. Comparison of different reconstitution methods after evaporation of 0.5 mL.





## GC Analysis of Pesticide Residues in Hops on SLB®-5ms after Extraction and Cleanup using Supel™ QuE Verde

There are a variety of pesticides used in the commercial production of hops. Hops contain essential oils, acids, pigments, and other compounds that present potential interferences. A standard approach to cleanup of hop extracts is the use of a PSA/C18/GCB sorbent blend. PSA removes acidic interferences while C18 removes hydrophobic interferences. GCB removes pigment. The compounds most effectively removed by commonly used GCB have planar structures. Any targeted pesticides with planar structures can also be retained.

In this application, a new sorbent blend, Supel QuE Verde, is used for QuEChERS cleanup in the analysis of pesticide residues in hops. Supel QuE Verde consists of PSA, Z-Sep+ and an improved GCB, which provides improved/superior recoveries of planar pesticides compared to traditional while still maintaining sufficient removal of the chlorophyll background.

<b>Sample/Matrix</b>	Ground hop pellets, 1 gm, spiked at 50 ng/g
<b>Sample Preparation:</b>	
<b>Extraction Process</b>	1g of ground sample. Add 10 mL of water to sample. Allow to sit for 30 minutes. Add 10 mL of acetonitrile containing 1% acetic acid; (shake for 10 min at 2500 rpm on automated shaker; Add contents of Supel QuE Citrate extraction tube (55227-U); shake for 1 minute; centrifuge at 5000 rpm for 5 min; remove supernatant)
<b>Cleanup Process</b>	Transfer 1 mL of supernatant to the Supel QuE Verde cleanup tube, 2 mL (55447-U); shake for 2 min; centrifuge at 5000 rpm for 5 min; remove supernatant; Place supernatant in amber low adsorption vial for GC/MS/MS analysis

### Results

Compound	% Rec	% RSD
Dichlorvos	68%	12%
Propoxur	86%	14%
Diazinon	101%	10%
Metalaxyl	83%	16%
Malathion	91%	13%
Chloropyrifos	97%	13%
Fipronil	89%	15%


Compound	% Rec	% RSD
Tetrachlorvinphos	82%	12%
Paclobutrazol	96%	10%
Myclobutanil	79%	7%
Pymetrozine	78%	7%
Trifloxystrobin	93%	9%
Quinoxifen	64%	23%
Spiromesifens	91%	19%

Compound	% Rec	% RSD
Piperonyl butoxide	90%	10%
Methoxychlor	101%	18%
Permethrins	76%	8%
Boscalid	76%	24%
Cypermethrins	87%	14%
Pyraclostrobin	107%	17%
Dimethomorph	58%	15%

### Analysis of a Standard Mix

column SLB-5ms, 30 m × 0.25 mm I.D. × 0.25 μm (28471-U)  
 oven 50 °C (2 min), 8 °C/min to 320 °C (5 min)  
 inj. temp. 250 °C  
 detector MS/MS  
 MSD interface 320 °C  
 carrier gas helium, 1.4 mL/min constant flow  
 injection 1 μL, splitless (splitter open at 0.75 min)  
 liner 4 mm I.D. FocusLiner with taper

- Dichlorvos
- Propoxur
- Diazinon
- Metalaxyl
- Propargite
- Malathion
- Chloropyrifos
- Fipronil
- Captan
- Tetrachlorvinphos
- Paclobutrazol
- Imazalil
- Pymetrozine
- Myclobutanil
- Trifloxystrobin
- Quinoxifen
- Tebuconazole
- Resmethrin isomers
- Piperonyl butoxide
- Spiromesifen isomers
- Tetramethrin isomers
- Phosmet
- Bifenthrin
- Fenoxycarb
- Etoxazole
- Methoxychlor
- γ-Cyhalothrin
- Spirodiclofen
- Permethrin isomers
- Cyfluthrin isomers
- Boscalid
- Dimethomorph isomers



## Sigma-Aldrich® Worldwide Offices

### Argentina

Free Tel: 0810 888 7446  
Tel: (+54) 11 4556 1472  
Fax: (+54) 11 4552 1698

### Australia

Free Tel: 1800 800 097  
Free Fax: 1800 800 096  
Tel: (+61) 2 9841 0555  
Fax: (+61) 2 9841 0500

### Austria

Tel: (+43) 1 605 81 10  
Fax: (+43) 1 605 81 20

### Belgium

Tel: (+32) 3 899 13 01  
Fax: (+32) 3 899 13 11

### Brazil

Free Tel: 0800 701 7425  
Tel: (+55) 11 3732 3100  
Fax: (+55) 11 5522 9895

### Canada

Free Tel: 1800 565 1400  
Free Fax: 1800 265 3858  
Tel: (+1) 905 829 9500  
Fax: (+1) 905 829 9292

### Chile

Tel: (+56) 2 495 7395  
Fax: (+56) 2 495 7396

### People's Republic of China

Free Tel: 800 819 3336  
Tel: (+86) 21 6141 5566  
Fax: (+86) 21 6141 5567

### Czech Republic

Tel: (+420) 246 003 200  
Fax: (+420) 246 003 291

### Denmark

Tel: (+45) 43 56 59 00  
Fax: (+45) 43 56 59 05

### Finland

Tel: (+358) 9 350 9250  
Fax: (+358) 9 350 92555

### France

Free Tel: 0800 211 408  
Free Fax: 0800 031 052  
Tel: (+33) 474 82 28 88  
Fax: (+33) 474 95 68 08

### Germany

Free Tel: 0800 51 55 000  
Free Fax: 0800 64 90 000  
Tel: (+49) 89 6513 0  
Fax: (+49) 89 6513 1169

### Hungary

Tel: (+36) 1 235 9055  
Fax: (+36) 1 235 9068

### India

#### Telephone

Bangalore: (+91) 80 6621 9400  
New Delhi: (+91) 11 4358 8000  
Mumbai: (+91) 22 4087 2364  
Pune: (+91) 20 4146 4700  
Hyderabad: (+91) 40 3067 7450  
Kolkata: (+91) 33 4013 8000

#### Fax

Bangalore: (+91) 80 6621 9550  
New Delhi: (+91) 11 4358 8001  
Mumbai: (+91) 22 2579 7589  
Pune: (+91) 20 4146 4777  
Hyderabad: (+91) 40 3067 7451  
Kolkata: (+91) 33 4013 8016

### Ireland

Free Tel: 1800 200 888  
Free Fax: 1800 600 222  
Tel: +353 (0) 402 20370  
Fax: + 353 (0) 402 20375

### Israel

Free Tel: 1 800 70 2222  
Tel: (+972) 8 948 4222  
Fax: (+972) 8 948 4200

### Italy

Free Tel: 800 827 018  
Tel: (+39) 02 3341 7310  
Fax: (+39) 02 3801 0737

### Japan

Tel: (+81) 3 5796 7300  
Fax: (+81) 3 5796 7315

### Korea

Free Tel: (+82) 80 023 7111  
Free Fax: (+82) 80 023 8111  
Tel: (+82) 31 329 9000  
Fax: (+82) 31 329 9090

### Luxembourg

Tel: (+32) 3 899 1301  
Fax: (+32) 3 899 1311

### Malaysia

Tel: (+60) 3 5635 3321  
Fax: (+60) 3 5635 4116

### Mexico

Free Tel: 01 800 007 5300  
Free Fax: 01 800 712 9920  
Tel: (+52) 722 276 1600  
Fax: (+52) 722 276 1601

### The Netherlands

Tel: (+31) 78 620 5411  
Fax: (+31) 78 620 5421

### New Zealand

Free Tel: 0800 936 666  
Free Fax: 0800 937 777  
Tel: (+61) 2 9841 0555  
Fax: (+61) 2 9841 0500

### Norway

Tel: (+47) 23 17 60 00  
Fax: (+47) 23 17 60 10

### Poland

Tel: (+48) 61 829 01 00  
Fax: (+48) 61 829 01 20

### Portugal

Free Tel: 800 202 180  
Free Fax: 800 202 178  
Tel: (+351) 21 924 2555  
Fax: (+351) 21 924 2610

### Russia

Free Tel: 8 800 100 7425  
Tel: (+7) 495 621 5828  
Fax: (+7) 495 621 6037

### Singapore

Tel: (+65) 6779 1200  
Fax: (+65) 6779 1822

### Slovakia

Tel: (+421) 255 571 562  
Fax: (+421) 255 571 564

### South Africa

Free Tel: 0800 1100 75  
Free Fax: 0800 1100 79  
Tel: (+27) 11 979 1188  
Fax: (+27) 11 979 1119

### Spain

Free Tel: 900 101 376  
Free Fax: 900 102 028  
Tel: (+34) 91 661 99 77  
Fax: (+34) 91 661 96 42

### Sweden

Tel: (+46) 8 742 4200  
Fax: (+46) 8 742 4243

### Switzerland

Free Tel: 0800 80 00 80  
Free Fax: 0800 80 00 81  
Tel: (+41) 81 755 2511  
Fax: (+41) 81 756 5449

### Thailand

Tel: (+66) 2 126 8141  
Fax: (+66) 2 126 8080

### United Kingdom

Free Tel: 0800 717 181  
Free Fax: 0800 378 785  
Tel: (+44) 01747 833 000  
Fax: (+44) 01747 833 574

### United States

Toll-Free: 800 325 3010  
Toll-Free Fax: 800 325 5052  
Tel: (+1) 314 771 5765  
Fax: (+1) 314 771 5757

### Vietnam

Tel: (+84) 8 3516 2810  
Fax: (+84) 8 6258 4238

### Internet

[sigma-aldrich.com](http://sigma-aldrich.com)

Order/Customer Service: [sigma-aldrich.com/order](http://sigma-aldrich.com/order)  
Technical Service: [sigma-aldrich.com/techservice](http://sigma-aldrich.com/techservice)  
Development/Custom Manufacturing Inquiries **SAFC** [safcglobal@sial.com](mailto:safcglobal@sial.com)  
Safety-related Information: [sigma-aldrich.com/safetycenter](http://sigma-aldrich.com/safetycenter)

3050 Spruce St.  
St. Louis, MO 63103  
(314) 771-5765  
[sigma-aldrich.com](http://sigma-aldrich.com)

©2016 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH, SAFC and SUPELCO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. Discovery, SLB, Ascentis and TraceCERT are registered trademarks of Sigma-Aldrich Co. LLC. PESTANAL is a registered trademark of Sigma-Aldrich Laborchemikalien GmbH. Supel, ENVI-Carb and Supelclean are trademarks of Sigma-Aldrich Co. LLC. Fused-Core is a registered trademark of Advanced Materials Technology, Inc. Benchmark is a trademark of Benchmark Scientific, Inc. Sigma-Aldrich, Supelco and SAFC brand products are sold by affiliated Sigma-Aldrich distributors. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at [www.sigmaaldrich.com](http://www.sigmaaldrich.com) and/or on the reverse side of the invoice or packing slip. Sigma-Aldrich Corp. is a subsidiary of Merck KGaA, Darmstadt, Germany.

RZO  
84444  
1026